

Food & beverage

VT306

VT306

Low Temperature Freezer –45°C

The VT low temperature freezers creates the possibility to maintain temperatures as low as -45°C. Supreme stability, reliability, user-friendliness and ease of cleaning make these freezers an ideal solution for laboratories and hospitals. R404a refrigerant, cyclopentane insulation combined with recyclable materials make the freezers extremely eco-friendly.

 **Energy Effecient**

 **Plug & Play**

 **Porthole**



VT306 | Gallery

VT306 |

Features	Value
No. of baskets	1 no.

VT306 |

Dimensions	Value
Dimension	885x1265x605
Weight	77/71

VT306 |

Cooling	Value
Refrigerant / amount	R290
Frequenze	50 Hz
Temperature range	-25 to -45 °C

